


CMR in Central vascular Imaging

Mohammed Chamsi Pasha, MD, FACC, FASE
 Assistant Professor, Weill Cornell Medical College
 Houston Methodist DeBakey Heart & Vascular Center

Methodist
 DeBakey Heart & Vascular Center



IG MEDICINE™

1


Contrast Enhanced MRA

- Gadolinium based contrast administered – first pass
- Imaging of entire 3D dataset can be done in 15-25 seconds (single breath hold)
- Advantages
 - Large FOV
 - High spatial resolution
 - High SNR

Methodist
 DeBakey Heart & Vascular Center

2

High resolution 3D angio



Methodist
 DeBakey Heart & Vascular Center

3


K-space Filling

- Ordering schemes
 - Linear
 - Centric
- Center : Image intensity
- Periphery : Fine detail

Methodist
 DeBakey Heart & Vascular Center

4

Location of Information in k-space



Center of k-space
 basic image contrast


Periphery of k-space
 edges, details

Methodist
 DeBakey Heart & Vascular Center

MRIquestions.com

5

Linear Ordering



Centric Ordering

Methodist
 DeBakey Heart & Vascular Center

6

Gadolinium contrast timing

- ⚡ Timing bolus
- ⚡ Fluoroscopic triggering
 - Automated
 - Manual
- ⚡ Time resolved

Methølist
DeBakey Heart
& Vascular Center

7

Timing Bolus

- ⚡ Test dose with 1-2 cc of contrast and 20-30 cc of saline flush
- ⚡ Measure transit time
- ⚡ Time full bolus such that middle part of infusion coincide with center of K space

Methølist
DeBakey Heart
& Vascular Center

8


Fluoroscopic Triggering

- ⚡ Begin contrast infusion
- ⚡ When contrast arrives to the vessel of interest, trigger 3D sequence

Methølist
DeBakey Heart
& Vascular Center

9


Care Bolus



- ⚡ Center of the k-space is measured as quickly as possible when the contrast agent reaches the region to be examined

Methølist
DeBakey Heart
& Vascular Center

10



Methølist
DeBakey Heart
& Vascular Center

JCMR 2012

11

Time Resolved MRA

- ⚡ Dynamic information about flow
- ⚡ Lower spatial resolution
- ⚡ Acquires multiple 3D data set every few seconds (5-8)
- ⚡ Outer part of K-space is sparsely sampled



Methølist
DeBakey Heart
& Vascular Center

12

VIBE - MRA


➤ **Volume inter-polated breath holding examination**

VIBE

MRI measurement technique: FLASH 3D imaging technique with reduced data acquisition time by using data interpolation, or partial Fourier techniques, or both, primarily for dynamic contrast-enhanced examinations of the abdomen.

Methølist
DeBakey Heart & Vascular Center

13




14

EKG gating

➤ **Motion is worst at the aortic root and along a mobile dissection flap, typically 250-300 msec**

Methølist
DeBakey Heart & Vascular Center

15



16


Visualization techniques

➤ **Maximal intensity projection**

Methølist
DeBakey Heart & Vascular Center

JCMR 2012

17



18

Comparison of Methods for Aorta Imaging

| Lesion | TTE | TOE | CT | MRI |
|--------------------------------|-----|-----|-----|-----|
| Ascending aortic dissection | ++ | +++ | +++ | +++ |
| Aortic arch dissection | + | + | +++ | +++ |
| Descending aortic dissection | + | +++ | +++ | +++ |
| Size | ++ | +++ | +++ | +++ |
| Mural thrombus | + | +++ | +++ | +++ |
| Intramural haematoma | + | +++ | ++ | +++ |
| Penetrating aortic ulcer | ++ | ++ | +++ | +++ |
| Involvement of aortic branches | + | (+) | +++ | +++ |

MethoList
DeBakey Heart
& Vascular Center

Erbel, European Heart Journal (2014) 35, 2873–2926

19

Choice of Imaging Modalities

Advantages of CT

- Universal availability
- Short scan times
- Higher spatial resolution
- Less cost
- Compatibility with ferromagnetic metals
- Coronary evaluation
- 2D/3D reconstructions
- Better for calcification


Advantages of CMR

- Tissue characterization, aortic wall imaging
- Dynamic imaging (dissection)
- Venous imaging
- AV pathology
- Pts with iodine allergy
- 4D flow
- Radiation free / young/pregnant/ACHD serial imaging

MethoList
DeBakey Heart
& Vascular Center

20

Standard Anatomic Landmarks of the Thoracoabdominal Aorta





| Anatomic Location |
|--|
| 1. Aortic sinuses of Valsalva |
| 2. Sinotubular junction |
| 3. Mid ascending aorta (midpoint in length between Nos. 2 and 4) |
| 4. Proximal aortic arch (aorta at the origin of the innominate artery) |
| 5. Mid aortic arch (between left common carotid and subclavian arteries) |
| 6. Proximal descending thoracic aorta (begins at the isthmus, approximately 2 cm distal to left subclavian artery) |
| 7. Mid descending aorta (midpoint in length between Nos. 6 and 8) |
| 8. Aorta at diaphragm (2 cm above the celiac axis origin) |
| 9. Abdominal aorta at the celiac axis origin |

Hiratzka. J Am Coll Cardiol. 2010;55(14):1509-1544

21

Effaced STJ Trileaflet AV



MethoList
DeBakey Heart
& Vascular Center

22

True short axis plane


- Diameter Is Measured Perpendicular to Axis of Blood Flow
- Double oblique method
- Outer to outer



MethoList
DeBakey Heart
& Vascular Center


23

MRA Descending aortic aneurysm Severe atherothrombosis




MethoList
DeBakey Heart
& Vascular Center

24



25




26

Coarctation Cine- CMR, non-contrast

- Discrete narrowing at the juxtaductal segment of the descending aorta
- Associations
 - HTN
 - Claudication
 - Turner's syndrome
 - ICH
 - BAV, Shon's complex

Methoist DeBakey Heart & Vascular Center

27




28

Markers of significant Coarac


- Reversal of intercostal flow
- Flow differential between proximal and distal end of the aorta

Methoist DeBakey Heart & Vascular Center

29



30



31

Role of Imaging in Aortic Dissection

- Determine extent of dissection
- Dimensions of aorta, branch vessels
- Viability of solid organs.

Life-threatening complications of dissection:

- Pericardial hemorrhage and tamponade
- Mediastinal hemorrhage
- Acute aortic valve insufficiency
- Coronary artery & carotid artery dissection
- End-organ malperfusion syndromes

Methoist
 DeBakey Heart & Vascular Center

32

Aortic Dissection Classification: 2 Systems

Classification of acute aortic syndrome

Stanford A
 Lesions affecting the ascending aorta with or without involvement of the descending aorta

DeBakey I
 Ascending and descending aorta


DeBakey II
 Ascending aorta only

Stanford B
 Lesions affecting the descending aorta only

DeBakey III
 Descending aorta only
 Left subclavian artery
 IIIa Affected region above aortic bifurcation
 IIIb Affected region extends below aortic bifurcation

Methoist
 DeBakey Heart & Vascular Center

33



34

CT/CMR Gold Standard for Diagnosis

Table 2. Reported Sensitivity and Specificity of Diagnostic Tools for Acute Aortic Syndrome^a

| Diagnostic Tool ^b | Studies, No. ^c | Patients, No. | Threshold | All AAS ^{d,e} | Sensitivity | Specificity |
|---|---------------------------|---------------|--------------------------------------|------------------------|----------------|-------------|
| CT ²⁷ | 1 | 49 | | 100 (86.3-100) | 100 | 100 |
| MRI ²⁹⁻³¹ | 3 | 116 | | 95.0-100 | 94.0-98.0 | |
| TTE ^{17,27,29,30,32,33} | 6 | 520 | | 86.0-100 | 90.0-100 | |
| TTE ^{34,35} | 2 | 228 | | 73.7-100 | 71.2-91.0 | |
| Intravascular ultrasound ³⁶ | 1 | 28 | | | | |
| D-dimer ³⁷⁻⁴² | 6 | 876 | >0.5-0.7 µg/mL | 51.7-100 | 32.8-89.2 | |
| Elastin degradation products ⁴³ | 1 | 609 | >3 SD above mean of healthy patients | 99.8 (99.1-100) | | |
| MMP 8/9 ⁴⁴ | 1 | 126 | >3.6 ng/mL | 100 (93.2-100) | 9.5 (3.9-18.5) | |
| Smooth muscle myosin heavy chain ⁴⁵ | 1 | 27 | >10 ng/mL | 90.0 (78.7-100) | 97.0 | |
| Soluble lectin-like oxidized LDLR 1 ⁴⁶ | 1 | 19 | >150 pg/mL | 89.5 | 94.3 | |

Methoist
 DeBakey Heart & Vascular Center

electrocardiographic gating is crucial for reducing motion artifact, particularly in the aortic root, ascending aorta, and proximal arch

35

Treatment Recommendations

Table 4. Treatment Recommendations for Acute Aortic Syndrome

| Stanford Type | ACCF/AHA/ATS/ACR/ASA/SCA/SCAI/SIR/STS/SVM 2010 ^{a,b} | Level (Grade) ^b | ESC 2014 ^{c,d} | Level (Grade) ^b | This Review | Level (Grade) ^b |
|----------------------------|---|---------------------------------|-------------------------|----------------------------|-------------------------|----------------------------|
| Aortic Dissection | | | | | | |
| A | Open surgical procedure | I (B) | Open surgical procedure | I (B) | Open surgical procedure | I (B) |
| B | Complicated ^e | Surgical procedure ^e | I (B) | TEVAR | I (C) | TEVAR |
| | Uncomplicated ^e | Medical therapy | I (B) | Medical therapy or TEVAR | I (C) or IIA (B) | Medical therapy or TEVAR |
| | | | | | | I (C) or IIA (C) |
| Intramural Hematoma | | | | | | |
| A | Complicated ^e | Open surgical procedure | IIA (C) | Open surgical procedure | I (C) | Open surgical procedure |
| | Uncomplicated ^e | Not mentioned | | Not mentioned | | Medical therapy |
| B | Complicated ^e | Surgical procedure ^e | IIA (C) | TEVAR | IIA (C) | TEVAR |
| | Uncomplicated ^e | Medical therapy | I (B) | Medical therapy | I (C) | Medical therapy |

Methoist
 DeBakey Heart & Vascular Center

Mussa. JAMA. 2016;316(7):754-763


36

Table 8. Suggested Follow-up of Aortic Pathologies After Repair or Treatment


| Pathology | Interval | Study |
|-----------------------|---|----------------------------------|
| Acute dissection | Before discharge, 1 mo, 6 mo, yearly | CT or MR, chest plus abdomen TTE |
| Chronic dissection | Before discharge, 1 y, 2 to 3 y | CT or MR, chest plus abdomen TTE |
| Aortic root repair | Before discharge, yearly | TTE |
| AVR plus ascending | Before discharge, yearly | TTE |
| Aortic arch | Before discharge, 1 y, 2 to 3 y | CT or MR, chest plus abdomen |
| Thoracic aortic stent | Before discharge, 1 mo, 2 mo, 6 mo, yearly Or 30 days* | CXR, CT, chest plus abdomen |
| Acute IMH/PAU | Before discharge, 1 mo, 3 mo, 6 mo, yearly | CT or MR, chest plus abdomen |

Hiratzka. J Am Coll Cardiol. 2010;55(14):1509-154


37




38




39



40




41



42

Aortic thrombus




The left image is a cross-sectional CT scan of the chest, showing a large, bright, filling defect within the aortic lumen, indicating a thrombus. The right image is a longitudinal CT scan of the aorta, showing a large, dark, filling defect within the lumen, also indicating a thrombus.

Methu
DeBakey Heart
& Vascular Center

43

T1 600 msec, Thrombus imaging



This T1-weighted MRI scan at 600 msec shows a cross-section of the heart. A large, bright, filling defect is visible in the left ventricle, which is characteristic of a thrombus. The surrounding myocardium and other cardiac structures are visible in various shades of gray.

Methodist
DeBakey Heart
& Vascular Center

44

Non contrast 3D SSFP


- Allows visualization of background tissues
- Does not require intravenous contrast agents
- Not dependent of accurate timing
- Useful for coronary imaging, and aorta

Method


DeBakey Heart
& Vascular Center

45

Coronary MRA



Resolution $1 \times 1 \times 1.5 \text{ mm}$
Acquisition time 10-15 minutes



T2 prep Nav F5 α -pulse SSFP data acquisition spoil

coronary stand-still

T2 prep Nav F5 α -pulse

Sheu SH, J MRI, 2003 10(5):977-902.


Metha
DeBakey & Vascular Center

46

47

Email:
machamsi-pasha@houstonmethodist.org

Twitter: @ChamsiPash



Methodist
DeBakey Heart
& Vascular Center

48


QISS (Quiescent-Interval single shot) Knees to toes

- ECG-triggered, single-shot 2D acquisition of one slice per heartbeat
- In-plane saturation to suppress background tissue and a tracking saturation pulse to suppress venous signal prior to a quiescent inflow period.

Methoist
DeBakey Heart
& Vascular Center

49

How it works



Methoist
DeBakey Heart
& Vascular Center

50



Methoist
DeBakey Heart
& Vascular Center

51